

Phyhd1 Rabbit Pab For Wb, If/Icc, Elisa - Q5Sre7

Item Number: CM0026530



Introduction

Rabbit polyclonal antibody against human PHYHD1, a 2-oxoglutarate-dependent dioxygenase. Suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat. Ships on ice bag.

[Learn More](#)

Parameter	Value
Product Name	PHYHD1 Rabbit pAb
Remarks/Alias	PHYHD1
Species	Human
GeneID (Human)	254295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-291 of human PHYHD1 (NP_001094346.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5SRE7
Protein Weight	32kDa
Shipping	Ice bag