


Pi3 Kinase P85 Alpha Rabbit Pab - P27986

Item Number: CM0017124



Introduction

High-quality rabbit polyclonal antibody targeting human PI3 Kinase p85 alpha (PIK3R1) for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for PI3K signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 alpha Rabbit pAb
Remarks/Alias	p85; AGM7; GRB1; IMD36; p85alpha; p85-ALPHA; PI3 Kinase p85 alpha
Species	Human
Gene ID (Human)	5295
Gene ID	5295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-380 of human PI3 Kinase p85 alpha (NP_852664.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27986
Protein Weight	84kDa
Shipping	Ice bag