


Pias1 Rabbit Polyclonal Antibody For If/Icc And Elisa - 075925

Item Number: CM0028586



Introduction

Rabbit polyclonal antibody targeting human E3 SUMO-protein ligase PIAS1, validated for IF/ICC and ELISA. Ideal for studying SUMOylation and transcriptional regulation pathways.

[Learn More](#)

Parameter	Value
Product Name	PIAS1 Rabbit pAb
Remarks/Alias	GBP; ZMIZ3; DDXBP1; GU/RH-II; PIAS1
Species	Human
GeneID (Human)	8554
GeneID	8554
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 422-651 of human PIAS1 (NP_057250.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	075925
Protein Weight	72kDa
Shipping	Ice bag