

Pias2 Rabbit Monoclonal Antibody For Wb, Elisa - O75928

Item Number: CM0004247



Introduction

PIAS2 Rabbit monoclonal antibody validated for WB and ELISA, reacts with human, mouse, rat; targets E3 SUMO-protein ligase PIAS2 involved in transcriptional regulation. Ideal for studying STAT, p53, and steroid hormone signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PIAS2 Rabbit mAb
Remarks/Alias	DIP; MIZ1; SIZ2; ARIP3; PIASX; ZMIZ4; PIAS2
Species	Human
GeneID (Human)	9063
GeneID	9063
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human PIAS2 (NP_004662.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O75928
Protein Weight	45kDa/63kDa/68kDa
Shipping	Ice bag