

Pias3 Rabbit Pab For Wb, Ihc-P, Elisa - Q9Y6X2

Item Number: CM0030940



Introduction

Rabbit polyclonal antibody against human PIAS3, suitable for WB, IHC-P, ELISA. Detects endogenous PIAS3 in human, mouse, and rat samples. E3 SUMO ligase involved in transcriptional regulation and cancer.

[Learn More](#)

Parameter	Value
Product Name	PIAS3 Rabbit pAb
Remarks/Alias	ZMIZ5; PIAS3
Species	Human
GeneID (Human)	10401
GeneID	10401
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 452-628 of human PIAS3 (NP_006090.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6X2
Protein Weight	68kDa
Shipping	Ice bag