

Pigr Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - 070570

Item Number: CM0030305



Introduction

PIGR Rabbit pAb validated for WB, IF/ICC, IHC-P, ELISA. Detects mouse/rat polymeric immunoglobulin receptor involved in transcytosis and mucosal immunity, ~83 kDa. Ideal for immunodetection studies.

[Learn More](#)

Parameter	Value
Product Name	PIGR Rabbit pAb
Species	Mouse
GeneID (Human)	5284
GeneID	18703
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-645 of mouse PIGR (NP_035212.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	070570
Protein Weight	83 kDa
Shipping	Ice bag