

Pik3Cg Rabbit Polyclonal Antibody For Wb, Elisa - P48736

Item Number: CM0015395



Introduction

Rabbit polyclonal antibody targeting human PIK3CG (p110γ) for Western blot and ELISA. Recognizes the 126 kDa PI3Kγ catalytic subunit involved in GPCR signaling and immune cell migration.

[Learn More](#)

Parameter	Value
Product Name	PIK3CG Rabbit pAb
Remarks/Alias	PI3K, PIK3, IMD97, PI3CG, PI3Kgamma, p110gamma, p120-PI3K
Species	Human
Gene ID	5294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human PIK3CG (NP_002640.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P48736
Protein Weight	126kDa
Shipping	Ice bag