

Pin4 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Y237

Item Number: CM0002175



Introduction

Rabbit monoclonal antibody targeting human PIN4, validated for WB and ELISA. Detects endogenous levels of PIN4 across human, mouse, and rat. Ideal for nucleolar and chromatin biology studies.

[Learn More](#)

Parameter	Value
Product Name	PIN4 Rabbit mAb
Remarks/Alias	EPVH; PAR14; PAR17; hEPVH; hPar14; hPar17; PIN4
Species	Human
GeneID (Human)	5303
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 52-131 of human PIN4 (Q9Y237).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y237
Protein Weight	14kDa
Shipping	Ice bag