

Pip5K1A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99755

Item Number: CM0025349

Introduction

This PIP5K1A rabbit polyclonal antibody detects human and mouse PIP5K1A (63 kDa) in WB, IHC-P, and ELISA. PIP5K1A is a lipid kinase that produces PIP2, regulating signal transduction, actin dynamics, and cell migration. Ideal for cell signaling and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	PIP5K1A Rabbit pAb
Remarks/Alias	PIP5K1A
Species	Human
Gene ID (Human)	8394
Gene ID	8394
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 231-500 of human PIP5K1A (NP_001129109.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99755
Protein Weight	63kDa
Shipping	Ice bag