

Piwil2 Rabbit Pab - Piwil2

Item Number: CM0030004

Introduction

PIWIL2 Rabbit pAb is a rabbit polyclonal antibody raised against a recombinant fusion protein of mouse PIWIL2. Validated for WB and ELISA, it detects mouse and rat PIWIL2 protein, making it ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	PIWIL2 Rabbit pAb
Remarks/Alias	mili; Piwil1l; PIWIL2
Species	Mouse
Gene ID (Human)	55124
Gene ID (Mouse)	57746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 398-652 of mouse PIWIL2 (Q8CDG1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
Protein Weight	105kDa/109kDa
Shipping	Ice bag