

Piwil4 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q7Z3Z4

Item Number: CM0005630

Introduction

PIWIL4 Rabbit Monoclonal Antibody (HIWI2, MIWI2) validated for WB, IF/ICC, ELISA with human, mouse, rat reactivity. Targets 97kDa PIWIL4 involved in piRNA-mediated transposon silencing. Ideal for germline and oncology research.

[Learn More](#)

Parameter	Value
Product Name	PIWIL4 Rabbit mAb
Remarks/Alias	HIWI2; MIWI2; PIWIL4
Species	Human
GeneID (Human)	143689
GeneID	143689
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 272-383 of human PIWIL4 (NP_689644.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z3Z4
Protein Weight	97kDa
Shipping	Ice bag