

## Piwil4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q7Z3Z4

Item Number: CM0018573

### Introduction

Rabbit polyclonal antibody targeting human PIWIL4 (HIWI2/MIWI2). Suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Detects piRNA pathway protein involved in transposon silencing and spermatogenesis.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | PIWIL4 Rabbit pAb                                                                                                   |
| Remarks/Alias    | HIWI2; MIWI2; PIWIL4                                                                                                |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 143689                                                                                                              |
| GeneID           | 143689                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 260-460 of human PIWIL4 (NP_689644.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q7Z3Z4                                                                                                              |
| Protein Weight   | 97kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |