


Pkc Alpha Rabbit Mab For If/Icc, Elisa - P17252

Item Number: CM0007536



Introduction

PKC alpha Rabbit mAb is a monoclonal antibody specific for human protein kinase C alpha, validated for IF/ICC and ELISA. Ideal for studying cell proliferation, apoptosis, and cancer signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PKC alpha Rabbit mAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha
Species	Human
GeneID	5578
Immunogen	Recombinant protein of human PKC alpha
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag