

Pkc-Beta Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05771

Item Number: CM0029916



Introduction

High-quality rabbit polyclonal antibody targeting human PKC-beta (PRKCB) for WB, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat. Ideal for studying B-cell receptor signaling, apoptosis, and insulin signaling.

[Learn More](#)

Parameter	Value
Product Name	PKC-beta Rabbit pAb
Remarks/Alias	PKCB; PRKCB1; PRKCB2; PKCI(2); PKCbeta; PKC-beta
Species	Human
Gene ID	5579
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-340 of human PKC-beta (NP_002729.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05771
Protein Weight	77kDa
Shipping	Ice bag