


Pkc Delta Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q05655

Item Number: CM0013187



Introduction

Rabbit polyclonal antibody targeting PKC delta (nPKC-delta, Q05655) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, and rat. Suitable for apoptosis, cancer, and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	PKC delta Rabbit pAb
Remarks/Alias	MAY1; PKCD; ALP53; CVID9; nPKC-delta; PKC delta
Species	Human
Gene ID (Human)	5580
Gene ID	5580
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human PKC delta (NP_006245.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05655
Protein Weight	78kDa
Shipping	Ice bag