


Pklr Rabbit Polyclonal Antibody For Wb And Elisa - P30613

Item Number: CM0016282



Introduction

PKLR rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human, mouse, and rat PKLR. Suitable for glycolysis and cancer metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PKLR Rabbit pAb
Remarks/Alias	PK1; PKL; RPK; PKRL; PKLR
Species	Human
GeneID (Human)	5313
GeneID	5313
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 446-556 of human PKLR (NP_000289.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30613
Protein Weight	62kDa
Shipping	Ice bag