


Pkm1-Specific Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P14618

Item Number: CM0002518



Introduction

PKM1-specific rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Targets the adult isoform of pyruvate kinase, involved in glycolysis. Reacts with human, mouse, rat. Suitable for metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PKM1-specific Rabbit mAb
Remarks/Alias	PK3; TCB; p58; OIP3; PKM2; CTHBP; THBP1; HEL-S-30; PKM1-specific
Species	Human
GeneID (Human)	5315
GeneID	5315
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PKM1 (NP_001193728.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14618
Protein Weight	58 kDa
Shipping	Ice bag