


Pkm1-Specific Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14618

Item Number: CM0022070



Introduction

Rabbit polyclonal antibody specifically recognizing human PKM1 (pyruvate kinase M1 isoform). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for studying glycolysis and energy metabolism in adult tissues.

[Learn More](#)

Parameter	Value
Product Name	PKM1-specific Rabbit pAb
Remarks/Alias	PK3; TCB; p58; OIP3; PKM2; CTHBP; THBP1; HEL-S-30; PKM1-specific
Species	Human
GeneID	5315
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PKM1-specific (NP_001193728.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14618
Protein Weight	58kDa
Shipping	Ice bag