


Plod3 Rabbit Mab - O60568

Item Number: CM0002894



Introduction

High-specificity rabbit monoclonal antibody targeting human PLOD3/LH3, validated for WB, IHC-P, and ELISA. Ideal for studying collagen lysine hydroxylation, glycosylation, and basement membrane biology.

[Learn More](#)

Parameter	Value
Product Name	PLOD3 Rabbit mAb
Remarks/Alias	LH3; PLOD3
Species	Human
Gene ID (Human)	8985
Gene ID	8985
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 25-306 of human PLOD3 (NP_001075.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O60568
Protein Weight	85kDa
Shipping	Ice bag