

Plxis Motif Rabbit Pab

Item Number: CM0020705

Introduction

Rabbit polyclonal antibody targeting the plxIS motif, derived from a synthetic peptide. Validated for Western blot (WB) and ELISA applications with human cross-reactivity. Ideal for research on proteins containing the plxIS motif.

[Learn More](#)

Parameter	Value
Product Name	plxIS motif Rabbit pAb
Immunogen	A synthesized peptide derived from plxIS motif
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
Shipping	Ice bag