

Pmpcb Rabbit Polyclonal Antibody For Western Blot And Elisa - 075439

Item Number: CM0026606

Introduction

Rabbit polyclonal antibody against mitochondrial-processing peptidase subunit beta (PMPCB) for Western blot and ELISA applications. Recognizes human, mouse, and rat PMPCB; suitable for mitochondrial protein import and processing research.

[Learn More](#)

Parameter	Value
Product Name	PMPCB Rabbit pAb
Remarks/Alias	MAS1, MPPB, P-52, MPP11, MPPP52, Beta-MPP, PMPCB
Species	Human
Gene ID (Human)	9512
Gene ID	9512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-380 of human PMPCB (NP_004270.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075439
Protein Weight	54kDa
Shipping	Ice bag