


Pms2 Rabbit Monoclonal Antibody For Wb, Elisa - P54278

Item Number: CM0010028



Introduction

High-specificity rabbit monoclonal antibody targeting human PMS2, validated in Western blot and ELISA. Essential for DNA mismatch repair studies and Lynch syndrome research. Detects a 96 kDa band (UniProt P54278).

[Learn More](#)

Parameter	Value
Product Name	PMS2 Rabbit mAb
Remarks/Alias	MLH4, PMS-2, PMSL2, HNPCC4, LYNCH4, MMRC54, PMS2CL, PMS2
Species	Human
GeneID (Human)	5395
GeneID	5395
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PMS2 (P54278)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P54278
Protein Weight	96kDa
Shipping	Ice bag