


Pnmt Rabbit Pab - P11086

Item Number: CM0030876



Introduction

PNMT rabbit polyclonal antibody for WB and ELISA applications. Detects human and mouse phenylethanolamine N-methyltransferase (31 kDa), a key enzyme in catecholamine biosynthesis. Suitable for IVD research.

[Learn More](#)

Parameter	Value
Product Name	PNMT Rabbit pAb
Remarks/Alias	PENT, PNMTase, PNMT
Species	Human
GeneID (Human)	5409
GeneID	5409
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-282 of human PNMT (NP_002677.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11086
Protein Weight	31kDa
Shipping	Ice bag