

Pnn Rabbit Pab For Wb, Elisa - Q9H307

Item Number: CM0022751

Introduction

PNN Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human and rat PNN, a transcriptional activator and splicing factor involved in cell adhesion and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	PNN Rabbit pAb
Remarks/Alias	DRS; DRSP; SDK3; memA; PNN
Species	Human
GeneID (Human)	5411
GeneID	5411
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PNN (NP_002678.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H307
Protein Weight	82kDa
Shipping	Ice bag