

Pnpla3 Rabbit Polyclonal Antibody (Adiponutrin, Ipla2-Epsilon) For Wb And Elisa - Q9Nst1

Item Number: CM0027850

Introduction

Rabbit polyclonal antibody against human PNPLA3 (adiponutrin, iPLA2-epsilon) for Western blot and ELISA. Detects human, mouse, and rat PNPLA3; suitable for lipid metabolism and NAFLD research.

[Learn More](#)

Parameter	Value
Product Name	PNPLA3 Rabbit pAb
Remarks/Alias	ADPN; C22orf20; iPLA(2)epsilon; PNPLA3
Species	Human
Gene ID (Human)	80339
Gene ID	80339
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 377-481 of human PNPLA3 (NP_079501.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NST1
Protein Weight	53kDa
Shipping	Ice bag