


Polr2E Rabbit Pab - P19388

Item Number: CM0019385



Introduction

Rabbit polyclonal antibody against human POLR2E (RPB5/RPABC1), a common RNA polymerase subunit. Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	POLR2E Rabbit pAb
Remarks/Alias	RPB5; XAP4; RPABC1; hRPB25; hsRPB5; POLR2E
Species	Human
Gene ID (Human)	5434
Gene ID	5434
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human POLR2E (NP_002686.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19388
Protein Weight	25kDa
Shipping	Ice bag