


Polr2l Rabbit Pab - P36954

Item Number: CM0020312



Introduction

DNA-directed RNA polymerase II subunit RPB9 (POLR2I) rabbit polyclonal antibody for Western blot and ELISA. Cross-reacts with human and mouse; integrates into Pol II transcription research.

[Learn More](#)

Parameter	Value
Product Name	POLR2I Rabbit pAb
Remarks/Alias	RPB9; hRPB14.5; POLR2I
Species	Human
Gene ID (Human)	5438
Gene ID	5438
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human POLR2I (NP_006224.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P36954
Protein Weight	15kDa
Shipping	Ice bag