


Pp2A Alpha/Beta Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - P67775 / P62714

Item Number: CM0009734



Introduction

PP2A alpha/beta rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat PP2A catalytic subunits, essential for serine/threonine phosphatase activity and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	PP2A alpha/beta Rabbit mAb
Remarks/Alias	PP2Abeta; PP2CB; PP2A alpha + beta
Species	Human
GeneID (human)	5515/5516
GeneID	5515/5516
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 200-309 of human PP2A alpha + beta (P62714)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67775/P62714
Protein Weight	36kDa
Shipping	Ice bag