

Ppa2 Rabbit Pab For Wb, If/Icc, Elisa - Q9H2U2

Item Number: CM0018826

Introduction

PPA2 Rabbit Polyclonal Antibody validated for Western Blot, Immunofluorescence, and ELISA. Cross-reacts with human, mouse, rat PPA2. Essential tool for studying mitochondrial pyrophosphatase and cellular energy metabolism.

[Learn More](#)

Parameter	Value
Product Name	PPA2 Rabbit pAb
Remarks/Alias	SCFI; SCFAI; HSPC124; SID6-306; PPA2
Species	Human
GeneID (Human)	27068
GeneID	27068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-100 of human PPA2 (NP_008834.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2U2
Protein Weight	38kDa
Shipping	Ice bag