

Ppar α Rabbit Pab - Q07869

Item Number: CM0028140



Introduction

PPAR α Rabbit pAb validated for WB, ELISA. Recognizes human, mouse, rat PPAR α . Ideal for lipid metabolism and nuclear receptor research. Reliable IVD raw material from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	PPAR α Rabbit pAb
Remarks/Alias	PPAR; NR1C1; hPPAR; PPARalpha; PPAR-alpha; PPAR α
Species	Human
GeneID (Human)	5465
GeneID	5465
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human PPAR α (NP_005027.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q07869
Protein Weight	52kDa
Shipping	Ice bag