

Ppbp Rabbit Mab For Wb, If/Icc, Elisa - P02775

Item Number: CM0006050



Introduction

High-quality PPBP Rabbit Monoclonal Antibody from CamelBio, validated for WB, IF/ICC, ELISA. Targets human platelet basic protein (CXCL7/NAP-2). Ideal for IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	PPBP Rabbit mAb
Remarks/Alias	PBP; TC1; TC2; TGB; LDGF; MDGF; TGB1; B-TG1; CTAP3; CXCL7; NAP-2; SCYB7; THBGB; LA-PF4; THBGB1; Beta-TG; CTAPIII; CTAP-III; PPBP
Species	Human
GeneID (human)	5473
GeneID	5473
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-128 of human PPBP (NP_002695.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P02775
Protein Weight	14kDa
Shipping	Ice bag