


Ppih Rabbit Pab - O43447

Item Number: CM0018600



Introduction

Recombinant rabbit polyclonal antibody targeting human PPIH (CypH), validated for WB and ELISA. Cross-reacts with mouse and rat. PPIH is a peptidyl-prolyl isomerase and spliceosome component, ideal for splicing and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	PPIH Rabbit pAb
Remarks/Alias	CYPH; CYP-20; USA-CYP; SnuCyp-20; PPIH
Species	Human
Gene ID (Human)	10465
Gene ID	10465
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43447
Protein Weight	19kDa
Shipping	Ice bag