


Ppil1 Rabbit Pab - Q9Y3C6

Item Number: CM0015450



Introduction

PPIL1 Rabbit pAb is a polyclonal antibody validated for Western blot and ELISA. It recognizes human PPIL1 (CYPL1), a spliceosome-associated peptidyl-prolyl isomerase involved in mRNA splicing and brain development. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	PPIL1 Rabbit pAb
Remarks/Alias	CYPL1; PCH14; hCyPX; PPlase; CGI-124; PPIL1
Species	Human
GeneID (Human)	51645
GeneID	51645
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-166 of human PPIL1 (NP_057143.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3C6
Protein Weight	18kDa
Shipping	Ice bag