


Ppil2 Rabbit Pab For Wb, Elisa - Q13356

Item Number: CM0027240



Introduction

High-quality PPIL2 rabbit polyclonal antibody validated for Western blot and ELISA applications, with cross-reactivity to human and mouse. Ideal for studying ubiquitin ligase function in splicing and protein degradation.

[Learn More](#)

Parameter	Value
Product Name	PPIL2 Rabbit pAb
Remarks/Alias	CYC4; CYP60; UBOX7; Cyp-60; hCyp-60; PPIL2
Species	Human
GeneID (Human)	23759
GeneID	23759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PPIL2 (NP_680481.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13356
Protein Weight	59kDa
Shipping	Ice bag