


Ppil4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wua2

Item Number: CM0024567



Introduction

PPIL4 rabbit polyclonal antibody for detection of the 57kDa Peptidyl-prolyl cis-trans isomerase-like 4 in human samples. Validated for Western blot (WB), immunofluorescence (IF/ICC), and ELISA.

[Learn More](#)

Parameter	Value
Product Name	PPIL4 Rabbit pAb
Remarks/Alias	HDCME13P; PPIL4
Species	Human
GeneID (Human)	85313
GeneID	85313
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-492 of human PPIL4 (NP_624311.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8WUA2
Protein Weight	57kDa
Shipping	Ice bag