

Ppm1G Rabbit Mab For Wb And Elisa - O15355

Item Number: CM0003623

Introduction

High-quality rabbit monoclonal antibody targeting human PPM1G (PP2Cγ). Validated for WB and ELISA, this antibody detects a 59kDa phosphatase involved in cellular signaling. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	PPM1G Rabbit mAb
Remarks/Alias	PP2CG; PPP2CG; PP2CGAMMA; PPM1G
Species	Human
GeneID (Human)	5496
GeneID	5496
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-450 of human PPM1G (O15355)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15355
Protein Weight	59kDa
Shipping	Ice bag