

Prmt1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99873

Item Number: CM0017484



Introduction

PRMT1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Targets Protein arginine N-methyltransferase 1 (UniProt: Q99873). Ideal for epigenetic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PRMT1 Rabbit pAb
Remarks/Alias	ANM1; HCPI1; IR1B4; HRMT1L2; PRMT1
Species	Human
GeneID (Human)	3276
GeneID	3276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-371 of human PRMT1 (NP_001527.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99873
Protein Weight	42kDa
Shipping	Ice bag