


Pro-Glucagon Rabbit Pab For Wb, If-P, Ihc-P, Elisa - P01275

Item Number: CM0015629



Introduction

Rabbit polyclonal antibody against Pro-glucagon, the precursor of glucagon, GLP-1, and GLP-2. Applications include WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for glucose metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Pro-glucagon Rabbit pAb
Remarks/Alias	GLP1; GLP2; GRPP; GLP-1; Pro-glucagon
Species	Human
Gene ID	2641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-89 of human Pro-glucagon (NP_002045.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01275
Protein Weight	21kDa
Shipping	Ice bag