

Prodh Rabbit Pab For Wb And Elisa - 043272

Item Number: CM0028298

Introduction

Polyclonal rabbit anti-PRODH antibody for Western blot and ELISA. Detects endogenous proline dehydrogenase 1 from human, mouse, and rat. The target is a mitochondrial enzyme central to proline metabolism and associated with cancer and schizophrenia.

[Learn More](#)

Parameter	Value
Product Name	PRODH Rabbit pAb
Remarks/Alias	POX; PIG6; HSPOX2; PRODH1; PRODH2; TP53I6; PRODH
Species	Human
Gene ID (Human)	5625
Gene ID	5625
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 321-600 of human PRODH (NP_057419.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43272
Protein Weight	68kDa
Shipping	Ice bag