


Proteasome 20S Lmp7 Rabbit Pab - P28062

Item Number: CM0021742



Introduction

Rabbit polyclonal antibody against human Proteasome 20S LMP7 subunit (PSMB8, 30kDa). Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for immunoproteasome and antigen processing research.

[Learn More](#)

Parameter	Value
Product Name	Proteasome 20S LMP7 Rabbit pAb
Remarks/Alias	JMP; ALDD; LMP7; NKJO; D6S216; PRAAS1; PSMB5i; RING10; D6S216E; Proteasome 20S LMP7
Species	Human
GeneID (human)	5696
GeneID	5696
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 167-276 of human Proteasome 20S LMP7. (NP_683720.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28062
Protein Weight	30kDa
Shipping	Ice bag