

Proz Rabbit Pab - P22891

Item Number: CM0023966



Introduction

Polyclonal antibody targeting vitamin K-dependent protein Z. Validated for Western blot and ELISA; cross-reacts with human, mouse, and rat. Ideal for coagulation and hemostasis studies.

[Learn More](#)

Parameter	Value
Product Name	PROZ Rabbit pAb
Remarks/Alias	PZ; PROZ
Species	Human
Gene ID (Human)	8858
Gene ID	8858
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 220-340 of human PROZ (NP_003882.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22891
Protein Weight	45kDa
Shipping	Ice bag