


Prpf31 Rabbit Mab For Wb, Elisa - Q8Wwy3

Item Number: CM0007011



Introduction

Recombinant rabbit monoclonal antibody targeting U4/U6 snRNP Prp31, a core spliceosome protein. Validated for WB and ELISA applications with human samples; 55 kDa.

[Learn More](#)

Parameter	Value
Product Name	PRPF31 Rabbit mAb
Remarks/Alias	RP11, PRP31, SNRNP61, NY-BR-99, PRPF31
Species	Human
GeneID (Human)	26121
GeneID	26121
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PRPF31 (Q8WWY3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WWY3
Protein Weight	55kDa
Shipping	Ice bag