


Prpf4 Rabbit Mab - O43172

Item Number: CM0006669



Introduction

PRPF4 Rabbit mAb is a monoclonal antibody targeting human PRPF4 (U4/U6 snRNP 60 kDa protein). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Suitable for splicing research and retinitis pigmentosa studies.

[Learn More](#)

Parameter	Value
Product Name	PRPF4 Rabbit mAb
Remarks/Alias	PRP4; RP70; HPRP4; Prp4p; HPRP4P; SNRNP60; PRPF4
Species	Human
Gene ID	9128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-522 of human PRPF4 (NP_004688.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43172
Protein Weight	58kDa
Shipping	Ice bag