


Prpf8 Rabbit Pab - Q6P2Q9

Item Number: CM0025533



Introduction

Rabbit polyclonal antibody targeting human PRPF8, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for spliceosome research.

[Learn More](#)

Parameter	Value
Product Name	PRPF8 Rabbit pAb
Remarks/Alias	PRP8; RP13; HPRP8; PRPC8; SNRNP220; PRPF8
Species	Human
GeneID (human)	10594
GeneID	10594
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PRPF8 (NP_006436.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P2Q9
Protein Weight	274kDa
Shipping	Ice bag