

Prrt2 Rabbit Pab For Wb And Elisa - Q7Z6L0

Item Number: CM0024422

Introduction

PRRT2 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human and mouse PRRT2, a regulator of neurotransmitter release involved in synaptic transmission. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PRRT2 Rabbit pAb
Remarks/Alias	PKC; EKD1; ICCA; BFIC2; BFIS2; DSPB3; DYT10; FICCA; IFITMD1; PRRT2
Species	Human
GeneID (human)	112476
GeneID	112476
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human PRRT2 (NP_660282.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z6L0
Protein Weight	35kDa
Shipping	Ice bag