

Prss2 Rabbit Polyclonal Antibody For Wb And Elisa - P07478

Item Number: CM0027671



Introduction

PRSS2 rabbit polyclonal antibody for Western blot and ELISA, reactive with human and mouse. PRSS2 (Trypsin-2) is a secreted serine protease involved in defensin processing in the ileum. Ideal for digestive enzyme and Paneth cell research.

[Learn More](#)

Parameter	Specification
Product Name	PRSS2 Rabbit pAb
Remarks/Alias	TRY2; TRY8; TRYP2; PRSS2
Species	Human
Gene ID (Human)	5645
Gene ID	5645
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-247 of human PRSS2 (NP_002761.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07478
Protein Weight	26kDa
Shipping	Ice bag