

Pseudouridine / 5-Ribosyluracil Rabbit Monoclonal Antibody For Db And Elisa Applications

Item Number: CM0000045

Introduction

Specific rabbit monoclonal antibody against pseudouridine (5-ribosyluracil) for dot blot and ELISA applications. Ideal for epitranscriptomics research and RNA modification analysis. Available from CamelBio as a reliable IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	Pseudouridine / 5-ribosyluracil Rabbit mAb
Gene ID (CAS)	CAS: 1445-07-4
Immunogen	Chemical compounds corresponding to Pseudouridine / 5-ribosyluracil / Y
Source	Rabbit
Category	Other Modified Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag