


Psma7 Rabbit Pab - O14818

Item Number: CM0025814



Introduction

PSMA7 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets proteasome subunit alpha type-7, 28kDa, involved in protein degradation and cellular regulation. Ideal for proteasome research.

[Learn More](#)

Parameter	Specification
Product Name	PSMA7 Rabbit pAb
Remarks/Alias	C6; HSPC; RC6-1; XAPC7; HEL-S-276; PSMA7
Species	Human
GeneID (Human)	5688
GeneID	5688
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-248 of human PSMA7 (NP_002783.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14818
Protein Weight	28kDa
Shipping	Ice bag