

Psmc3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17980

Item Number: CM0020199



Introduction

Rabbit anti-PSMC3 polyclonal antibody validated for Western blot, immunohistochemistry, and ELISA. Detects human, mouse, and rat PSMC3, a key 26S proteasome regulatory subunit. Ideal for ubiquitin-proteasome pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PSMC3 Rabbit pAb
Remarks/Alias	RPT5; TBP1; DCIDP; PSMC3
Species	Human
Gene ID (Human)	5702
Gene ID	5702
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human PSMC3 (NP_002799.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17980
Protein Weight	49kDa
Shipping	Ice bag