


Psmd12 Rabbit Polyclonal Antibody For Wb, Elisa - 000232

Item Number: CM0029241



Introduction

PSMD12 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat PSMD12 (26S proteasome non-ATPase regulatory subunit 12), a key component of the ubiquitin-proteasome system. Ideal for proteasome research.

[Learn More](#)

Parameter	Value
Product Name	PSMD12 Rabbit pAb
Remarks/Alias	p55; Rpn5; STISS; PSMD12
Species	Human
Gene ID (Human)	5718
Gene ID	5718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 197-456 of human PSMD12 (NP_002807.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00232
Protein Weight	53kDa
Shipping	Ice bag